

Hepatitis B e Antigen, Serum

Test ID: EAG

Explanation: Due to a method change, the following changes will be implemented on the effective date.

Current Method
Chemiluminescence Immunoassay

New Method
Electrochemiluminescence Immunoassay (ECLIA)

Current Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	28 days
	Ambient	24 hours
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	90 days
	Ambient	72 hours
	Refrigerated	6 days

Current Requested Minimum Volume
0.5 mL

New Requested Minimum Volume
0.6 mL

Questions

Contact Dunisha Messmer, Laboratory Resource Coordinator at 800-533-1710.