

## Cystinuria Gene Panel, Varies

**Test ID:** CYSGP

**Explanation:** On the effective date, this assay will no longer accept *any anticoagulant* for the whole blood specimen type and will be limited to the anticoagulants shown below. Additionally, the Report Available time will be updated.

| Current Specimen Required              |  |
|--|--|
| <b>Specimen Type:</b>                  | <b>Whole blood</b>   |
| <b>Preferred:</b>                      | Lavender top (EDTA) or yellow top (ACD)  |
| <b>Acceptable:</b>                     | <b>Any anticoagulant</b>   |
| <b>Specimen Volume:</b>                | 3 mL   |
| <b>Collection Instructions:</b>        | 1. Invert several times to mix blood<br>2. Send whole blood specimen in original tube.<br><b>Do not aliquot.</b> |
| <b>Specimen Stability Information:</b> | Ambient (preferred) 4 days/Refrigerated 14 days  |
| <b>Minimum Volume:</b>                 | 1 mL   |

| New Specimen Required                  |  |
|--|--|
| <b>Specimen Type:</b>                  | <b>Whole blood</b>   |
| <b>Required:</b>                       | Lavender top (EDTA) or yellow top (ACD)  |
| <b>Specimen Volume:</b>                | 3 mL   |
| <b>Collection Instructions:</b>        | 1. Invert several times to mix blood<br>2. Send whole blood specimen in original tube.<br><b>Do not aliquot.</b> |
| <b>Specimen Stability Information:</b> | Ambient (preferred) 4 days/Refrigerated 14 days  |
| <b>Minimum Volume:</b>                 | 1 mL   |

| Current Report Available |
|--------------------------|
| 28 to 42 days            |

| New Report Available |
|----------------------|
| 21-28 days           |

### Questions

Contact Michelle Rath, Laboratory Resource Coordinator at 800-533-1710.