

## Hepatitis B Virus Surface Antibody Screen, Qualitative/Quantitative, Serum

**Test ID:** HBBSN

**Explanation:** Due to manufacturer's modification of stability times, the stability information will be updated as indicated below.

Current Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	90 days
	Ambient	72 hours
	Refrigerated	6 days

New Specimen Stability		
Specimen	Temperature	Time
Serum SST	Frozen (preferred)	90 days
	Ambient	7 hours
	Refrigerated	6 days

### Questions

Contact Dunisha Messmer, Laboratory Resource Coordinator at 800-533-1710.