

Thyroid-stimulating Immunoglobulin (TSI)

Explanation: The following new orderable, referred to LabCorp Burlington, is available immediately.

Test ID: FFTSI

Methods:

IMMULITE 2000 TSI

Reference Values:

0.00-0.55 IU/L

Specimen Requirements:

Collection Container/Tube:

Preferred: Red-top

Acceptable: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 3 mL

Collection Instructions:

1. Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable.
2. Within 2 hours of draw, centrifuge and aliquot 3 mL of serum into a plastic vial.
3. Send refrigerate.

Specimen Minimum Volume: 0.3 mL (Note: This volume does not allow for repeat testing.)

Specimen Stability Information:

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Ambient	14 days
	Frozen	14 days

Fee: \$85.00

CPT Code: 84445

Day(s) Setup: Monday through Friday

Report Available: 4 to 8 days

Questions:

Contact MCL Referrals Supervisor Amy Bluhm at 800-533-1710.