

Eculizumab, Serum

Test ID: ECULI

Explanation:

On the effective date, stability for ECULI will be updated as shown below.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	30 days
	Ambient	4 days
	Refrigerated	28 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Frozen (preferred)	28 days
	Ambient	28 days
	Refrigerated	28 days

Questions

Contact Joshua Yang, Laboratory Resource Coordinator at 800-533-1710.