

Sequential Maternal Screening, Part 2, Serum

Test ID: SEQB

Explanation:

On the aforementioned date, test Sequential Maternal Screening, Part 2, Serum will become obsolete.

Recommended Alternative Test:

Maternal Screening, Sequential, Specimen #2, Alpha Fetoprotein, hCG, Estriol, and Inhibin A

Test ID: ZW242 - 3000148

Methodology:

Quantitative Chemiluminescent Immunoassay

Specimen Requirements:

Preferred:	Serum Separator Tube (SST)
Acceptable:	Plain Red
Specimen Volume:	3 mL
Collection Instructions:	Separate from cells ASAP or within 2 hours of collection. Transfer 3 mL serum to an ARUP Standard Transport Tube.
Minimum Volume:	1 mL
Additional Information:	<p>Submit with Maternal Serum Testing History form. Specimen must be drawn between 14 weeks, 0 days and 24 weeks, 6 days gestation. The recommended time for maternal serum screening is 16 to 18 weeks gestation.</p> <p>Requires that a previous first trimester specimen, Maternal Screening, Sequential, Specimen #1, hCG, PAPP-A, NT (ARUP test code 3000146), has been performed.</p>

Specimen Stability Information:

Specimen Type	Temperature	Time
Serum	Refrigerate (preferred)	14 days
	Frozen	1 year
	Ambient	72 hours

CPT Code:

81511

Miscellaneous Referral Test Fee:

\$160.00

Questions

Contact Bethany Feind, Laboratory Resource Coordinator at 800-533-1710.