

Ethyl Glucuronide Screen, Random, Urine

Test ID: ETGX

Explanation: Clarifying language is being added to reference value field and ambient specimen stability has been extended.

Current Reference Value	New Reference Value
Negative Cutoff concentrations: 500 ng/mL	Negative
	Positive results are reported with a quantitative result.
	Cutoff Concentrations:
	Immunoassay screen: 500 ng/mL
	Liquid chromatography tandem mass spectrometry: Ethyl glucuronide: 250 ng/mL Ethyl sulfate: 100 ng/mL

Current Specimen Stability			
Specimen Type	Temperature	Time	
Urine	Refrigerated (preferred)	28 days	
	Ambient	72 hours	
	Frozen	28 days	

New Specimen Stability			
Specimen Type	Temperature	Time	
Urine	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Questions Contact Brandon DeBoom, Laboratory Resource Coordinator at 800-533-1710.