

Ethyl Glucuronide Screen with Reflex, Random, Urine

Test ID: ETGR

Explanation: Ambient specimen stability has been extended from 72 hours to 28 days.

Current Specimen Stability		
Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	28 days
	Ambient	72 hours
	Frozen	28 days

New Specimen Stability		
Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	28 days
	Ambient	28 days
	Frozen	28 days

Questions

Contact Brandon DeBoom, Laboratory Resource Coordinator at 800-533-1710.