

Methadone and Metabolites, Serum

Test ID: MDNS

Explanation: Clarifying language and cutoff concentrations are being added to the reference value field as shown below.

Current Reference Value	New Reference Value
Not established	Methadone: Therapeutic: 100-400 ng/mL Toxic: >2000 ng/mL
	EDDP: Not established
	Cutoff concentrations by liquid chromatography tandem mass spectrometry: Methadone: 10 ng/mL EDDP: 5.0 ng/mL

Questions

Contact Brandon DeBoom, Laboratory Resource Coordinator at 800-533-1710.