

Notification Date: May 20, 2025 Effective Date: May 28, 2025

# T-Cell Receptor Constant Beta Chain (TRBC1) Immunostain, Technical Component Only

Test ID: TRBC1

#### **Useful for:**

Determination of T-cell clonality in T-cell neoplasms

## **Reflex Tests:**

Test ID	Reporting Name	Available Separately	Always Performed
IHTOI	IHC Initial, Tech Only	No	No
IHTOA	IHC Additional, Tech Only	No	No

## Methods:

Immunohistochemistry (IHC)

## **Specimen Requirements:**

**Supplies:** Immunostain Technical Only Envelope (T693)

Specimen Type: Tissue

Container/Tube: Immunostain Technical Only Envelope (T693)

Preferred: Formalin-fixed, paraffin-embedded (FFPE) tissue block

OR

2 unstained positively charged glass slide (25 x 75 x 1mm) per test ordered;

sections 4-microns thick

## **Specimen Stability Information:**

Specimen Type	Temperature	Time
TECHONLY	Ambient (preferred)	
	Refrigerated	

## **Cautions:**

Age of a cut paraffin section can affect immunoreactivity. Stability thresholds vary widely among published literature and are antigen dependent. Best practice is for paraffin sections to be cut within 6 weeks.

The charge of glass slides can be affected by environmental factors and subsequently may alter slide staining. Sending unsuitable glass slides can result in inconsistent staining due to poor slide surface chemistry.

Best practices for storage of positively charged slides:

- -Minimize time slides are stored after being unpackaged
- -Limit exposure to high humidity and heat
- -Minimize exposure to plastics

#### **CPT Code:**

88342-TC, primary 88341-TC, if additional IHC

**Day(s) Performed:** Monday through Friday **Report Available:** 1 day (Digital Image)

#### **Questions**

Contact Richard Einerson/Sarah Turner, Laboratory Technologist Resource Coordinator at 800-533-1710.