

Iodine/Creatinine Ratio, Random, Urine

Test ID: IODCU

Explanation: Acceptable stability for IODCU specimens will change on the effective date.

Current Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated (preferred)	30 days
	Ambient	14 days
	Frozen	30 days

New Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated (preferred)	28 days
	Ambient	14 days
	Frozen	28 days

Questions

Contact Rebekah Walsh, Laboratory Resource Coordinator at 800-533-1710.