

Iodine, 24 Hour, Urine

Test ID: UIOD

Explanation: Acceptable stability for UIOD specimens will change on the effective date.

Current Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated (preferred)	146 days
	Ambient	146 days
	Frozen	146 days

New Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerated (preferred)	28 days
	Ambient	28 days
	Frozen	28 days

Questions

Contact Rebekah Walsh, Laboratory Resource Coordinator at 800-533-1710.