

Neonatal Alloimmune Thrombocytopenia (NAIT)-Serial Monitoring Maternal only

Explanation: The following new orderable, referred to Versiti Wisconsin, Inc., will be available on August 28, 2025.

Test ID: FNNAT

Methods:

Platelet Antibody Bead Array (PABA) and Flow Cytometry

Reference Values:

Negative

Specimen Requirements:

Specimen Type: Serum

Container/Tube: Red top tube

Specimen Volume: 5 mL

Collection Instructions:

1. Sample should be spun down and taken off the clot.
2. Send 5 mL serum. Ship Refrigerated.

Specimen Stability Information:

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	365 days

Fee:

\$ 314.00

CPT Code:

86022, 86023 x 2

Day(s) Setup: Varies

Report Available: 9-11 days

Questions:

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.