

## Boron, Serum/Plasma

**Test ID:** BOAC

**Explanation:** Per notification received from NMS Labs, effective October 6, 2025, Test ID BOAC, will experience the following changes.

Current Method
Inductively Coupled Plasma/Optical Emission Spectrometry (ICP/OES)

New Method
Inductively Coupled Plasma/Mass Spectrometry (ICP/MS)

Current Reference Value
Normally: Less than 100 mcg/L

New Reference Value
Typically less than 100 mcg/L

Current Analytic Time
9 days

New Analytic Time
4 to 9 days

### Questions:

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.