

Bath Salts Panel (Qualitative), Urine

Test ID: FBSQU

Explanation:

Effective November 3, 2025, Test ID FBSQU, referred to NMS Labs, will be obsolete and replaced with following test.

Recommended Alternative Test:

Miscellaneous NMS Testing

Test ID: ZW86

Referral Lab Code: 2626U

Test Name: NPS Stimulants/Bath Salts Panel, Urine

Methodology:

High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

Reference Values:

N,N-Dimethylpentylone: None detected

Pentylone: None detected

Ethylone: None detected

Butylone: None detected

N-ethyl Pentylone: None detected

Specimen Requirements:

Must submit one specimen per order. Specimens cannot be shared between multiple orders.

Specimen Type: Urine

Container/Tube: Plastic preservative-free urine container

Specimen Volume: 1 mL

Collection Instructions:

1. Collect 1 mL random urine without preservative.
2. Send specimen refrigerated.

Minimum Volume: 0.7 mL

Specimen Stability Information:

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	14 days
	Ambient	7 days
	Frozen	14 days

Fee: \$141.00

CPT Code: 80371

Day(s) Setup: Monday through Sunday

Report Available: 6 to 12 days

Questions

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.