

## 5'Nucleotidase

**Test ID:** F5NUL

**Explanation:**

Effective November 06, 2025, Test ID F5NUL, referred to ARUP Laboratories, will be obsolete due to availability of Mayo tests 5NTD.

**Recommended Alternative Test:**

## 5'-Nucleotidase, Serum

**Test ID:** 5NTD

**Methodology:**

Coupled Enzyme Kinetic Assay

**Reference Values:**

<4.6 U/L

For International System of Units (SI) conversion for Reference Values, see [www.mayocliniclabs.com/order-tests/si-unit-conversion.html](http://www.mayocliniclabs.com/order-tests/si-unit-conversion.html).

**Specimen Requirements:**

**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)

**Collection Container/Tube:**

**Preferred:** Serum gel

**Acceptable:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

**Collection Instructions:** Centrifuge and aliquot serum into plastic vial.

**Specimen Stability Information:**

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	14 days
	Frozen	90 days

**CPT Code:**

83915

**Day(s) Setup:** Monday through Saturday**Analytic Time:** 1 day**Questions**

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.