

Porphobilinogen, Quantitative, Random, Urine

Test ID: PBGU

Explanation: The reference value will be updated due to a conversion to alternative units of measure (mmol/mol Cr).

Current Reference Value
< or =1.3 mcmol/L

New Reference Value
< or =0.2 mmol/mol creatinine

Current Specimen Required
Patient Preparation: Patient must not consume any alcohol for at least 24 hours prior to collection. Supplies: Urine Container-Amber, 60 mL (T596) Specimen Volume: 20 mL Collection Instructions: 1. Collect a random urine specimen. 2. No preservative necessary but pH must be above 5.0. 3. Specimens should be protected from light and frozen immediately following collection

New Specimen Required
Patient Preparation: Patient must not consume any alcohol for at least 24 hours before specimen collection. Supplies: Urine Container-Amber, 60 mL (T596) Specimen Volume: 20 mL Collection Instructions: 1. Collect a random urine specimen. 2. No preservative necessary but pH must be above 5.0. If pH <5.0, specimen will be rejected. 3. Specimens should be protected from light and frozen immediately following collection.

Questions

Contact Melissa Tricker-Klar, Laboratory Resource Coordinator at 800-533-1710.