

## Pneumococcal Antibody Panel - PCV15 (15-serotype)

**Explanation:** The following new orderable, referred to Eurofins Viracor, will be available on January 19, 2026.

**Test ID:** FPA15

**Methods:**

Fluoroimmunoassay Multiplex Bead Assay

**Reference Values:**

Pneumococcal Antibody Panel - PCV15

Pneumo Ab Type 1: >1.3 ug/mL

Pneumo Ab Type 3: >1.3 ug/mL

Pneumo Ab Type 4: >1.3 ug/mL

Pneumo Ab Type 5: >1.3 ug/mL

Pneumo Ab Type 6A: >1.3 ug/mL

Pneumo Ab Type 6B: >1.3 ug/mL

Pneumo Ab Type 7F: >1.3 ug/mL

Pneumo Ab Type 9V: >1.3 ug/mL

Pneumo Ab Type 14: >1.3 ug/mL

Pneumo Ab Type 18C: >1.3 ug/mL

Pneumo Ab Type 19A: >1.3 ug/mL

Pneumo Ab Type 19F: >1.3 ug/mL

Pneumo Ab Type 22F: >1.3 ug/mL

Pneumo Ab Type 23F: >1.3 ug/mL

Pneumo Ab Type 33F: >1.3 ug/mL

**Specimen Requirements:**

**Collection Container/Tube:**

**Preferred:** Red top

**Acceptable:** Serum gel

**Submission Container/tube:** Plastic vial

**Specimen Volume:** 1 mL serum

**Collection Instructions:**

1. Centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send frozen.

**Specimen Minimum Volume**

Serum: 0.1 mL

**Specimen Stability Information:**

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	42 days
	Ambient	42 days
	Frozen	

**Fee:**  
\$140.00

**CPT Code:**  
86581

**Day(s) Setup:** Monday through Friday

**Report Available:** 9 to 11 days

**Questions**

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.