

Pneumococcal Antibody Panel (12 Serotype)

Test ID: FPNAP

Explanation:

Effective January 19, 2026, Test ID FPNAP, referred to Eurofins Viracor, will be obsolete and replaced with Test ID FPA15.

Recommended Alternative Test:

Pneumococcal Antibody Panel - PCV15 (15-serotype)

Test ID: FPA15

Methodology:

Fluoroimmunoassay Multiplex Bead Assay

Reference Values:

Pneumococcal Antibody Panel - PCV15

Pneumo Ab Type 1: >1.3 ug/mL

Pneumo Ab Type 3: >1.3 ug/mL

Pneumo Ab Type 4: >1.3 ug/mL

Pneumo Ab Type 5: >1.3 ug/mL

Pneumo Ab Type 6A: >1.3 ug/mL

Pneumo Ab Type 6B: >1.3 ug/mL

Pneumo Ab Type 7F: >1.3 ug/mL

Pneumo Ab Type 9V: >1.3 ug/mL

Pneumo Ab Type 14: >1.3 ug/mL

Pneumo Ab Type 18C: >1.3 ug/mL

Pneumo Ab Type 19A: >1.3 ug/mL

Pneumo Ab Type 19F: >1.3 ug/mL

Pneumo Ab Type 22F: >1.3 ug/mL

Pneumo Ab Type 23F: >1.3 ug/mL

Pneumo Ab Type 33F: >1.3 ug/mL

Specimen Requirements:**Collection Container/Tube:**

Preferred: Red top

Acceptable: Serum gel

Submission Container/tube: Plastic vial

Specimen Volume: 1 mL serum

Collection Instructions:

1. Centrifuge and aliquot 1 mL of serum into a plastic vial.
2. Send frozen.

Specimen Minimum Volume

Serum: 0.1 mL

Specimen Stability Information:

Specimen Type	Temperature	Time
Serum	Refrigerated (preferred)	42 days
	Ambient	42 days
	Frozen	

CPT Code:

86581

Day(s) Setup: Monday through Friday

Report Available: 9 to 11 days

Questions

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.