

## Lactoferrin, Fecal by ELISA

**Test ID:** FLACF

**Explanation:**

Effective January 15, 2026, Test ID FLACF, referred to ARUP Laboratories, will be obsolete due to availability of Mayo test LCTF.

**Recommended Alternative Test:**

## Lactoferrin, Feces

**Test ID:** LCTF

**Methodology:**

Enzyme-Linked Immunosorbent Assay (ELISA)

**Reference Values:**

< or =4.4 mcg/g (Normal)

>4.4 mcg/g (Abnormal)

Reference values apply to all ages.

**Specimen Requirements:**

**Supplies:** Stool container, Small (Random), 4 oz (T288)

**Container/Tube:** Stool container

**Specimen Volume:** 5 g

**Collection Instructions:**

1. Collect a fresh random fecal specimen, no preservatives.
2. If specimen is sent refrigerate, send immediately after collection.
3. If specimen cannot be sent immediately, freeze specimen and send frozen.

### Additional Information:

1. Separate specimens must be submitted when multiple tests are ordered, with the exception of ELASF / Pancreatic Elastase, Feces and CALPR / Calprotectin, Feces. If only a single specimen is collected, it must be split prior to transport.

2. Testing cannot be added on to a previously collected specimen.

**Specimen Minimum Volume:** 1 g

### Specimen Stability Information:

Specimen Type	Temperature	Time
Fecal	Frozen (preferred)	7 days
	Ambient	72 hours
	Refrigerated	72 hours

### CPT Code:

83631

**Day(s) Setup:** Tuesday

**Report Available:** 3 to 8 days

### Questions

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.