

## Thiothixene (Navane)

**Test ID:** FFTHI

**Explanation:**

Effective immediately, Test ID FFTHI, referred to Medtox Laboratories, will become obsolete due to discontinuation at the performing laboratory.

**Recommended Alternative Test:**

## Miscellaneous NMS Testing

**Test ID:** ZW86

**Referral Lab Code:** 4469SP

**Test Name:** Thiothixene, Serum/Plasma

**Methodology:**

High Performance Liquid Chromatography/ Tandem Mass Spectrometry (LC-MS/MS)

**Reference Values:**

Interpretive Report

**Specimen Requirements:**

\*\*\*Must submit one specimen per order. Specimens cannot be shared between multiple orders\*\*\*

Submit one of the following:

Collection Container/Tube: Red top tube

Submission Container/Tube: Plastic screw-capped vial

Specimen Volume: 1 mL Serum

Collection Instruction:

1. Promptly centrifuge and separate Serum into a plastic screw capped vial (preservative-free plastic container).
2. Ship Refrigerated

Collection Container/Tube: Lavender top tube (EDTA) or Pink top tube

Submission Container/Tube: Plastic screw-capped vial

Specimen Volume: 1 mL Plasma

Collection Instruction:

1. Promptly centrifuge and separate Plasma into a plastic screw capped vial (preservative-free plastic container).
2. Ship Refrigerated

Min Vol: 0.5 mL

### Specimen Stability Information:

Specimen Type	Temperature	Time
Varies	Refrigerated (preferred)	14 days
	Ambient	7 days
	Frozen	

### CPT Code:

80342

### Fee:

\$125.00

**Day(s) Setup:** Monday through Sunday

**Report Available:** 6 to 10 days

### Questions

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.