

Cortisol, Free, Random, Urine

Test ID: CRANR

Explanation: On the effective date, CRANR will become obsolete. Please order the recommended alternative test below.

Recommended Alternative Test:

Cortisol, Free, 24 Hour, Urine

Test ID: CORTU

Methodology:

Liquid Chromatography Tandem Mass Spectrometry (LC-MS/MS)

Reference Values:

0-2 years: Not established

3-8 years: 1.4-20 mcg/24 h

9-12 years: 2.6-37 mcg/24 h

13-17 years: 4.0-56 mcg/24 h

> or =18 years: 3.5-45 mcg/24 h

Specimen Requirements:

Container/Tube: Urine Tubes, 10-mL (T068)

Specimen Volume: 5 mL

Collection Instructions:

1. Collect urine for 24 hours.
2. Add 10 g of boric acid as preservative at start of collection.

Minimum Volume: 3 mL

Specimen Stability Information:

Specimen Type	Temperature	Time
Urine	Refrigerated (Preferred)	14 days

	Ambient	7 days
	Frozen	28 days

CPT Code:

82530

Day(s) Performed: Monday - Friday**Report Available:** 2 to 7 days**Questions**

Contact Joshua Yang, Laboratory Resource Coordinator at 800-533-1710.