

## Hypersensitivity Pneumonitis Panel

**Test ID:** FHSP

### Explanation:

Effective April 23, 2026, Test ID FHSP, referred to Eurofins Viracor Clinical Diagnostics, will be obsolete due to availability of Mayo test ID BHYS.

### Recommended Alternative Test:

## Hypersensitivity Pneumonitis Panel, Immunoglobulin G, Serum

**Test ID:** BHYS

### Methodology:

Fluorescence Enzyme Immunoassay (FEIA)

### Reference Values:

*Alternaria alternata*, IgG antibody: < or =19.0 mg/L  
*Aspergillus fumigatus*, IgG antibody: < or =102.0 mg/L  
*Aureobasidium pullulans*, IgG antibody: < or =16.0 mg/L  
*Laceyella sacchari*, IgG antibody: < or =45.0 mg/L  
*Micropolyspora faeni*, IgG antibody: < or =6.0 mg/L  
*Penicillium chrysogenum*, IgG antibody: < or =94.0 mg/L  
*Phoma betae*, IgG antibody: < or =16.0 mg/L  
*Trichoderma viride*, IgG antibody: < or =16.0 mg/L

### Specimen Requirements:

**Container/Tube:** Sarstedt Aliquot Tube, 5 mL (T914)

**Preferred:** Serum gel

**Acceptable:** Red top

**Specimen Volume:** 0.5 mL Serum

**Collection Information:** Centrifuge and aliquot serum into a plastic vial.

### Specimen Stability Information:

| <b>Specimen Type</b> | <b>Temperature</b>       | <b>Time</b> |
|----------------------|--------------------------|-------------|
| Serum                | Refrigerated (preferred) | 28 days     |
|                      | Frozen                   | 28 days     |

**CPT Code:**

86001 x 8

**Day(s) Setup:** Monday through Friday

**Report Available:** 2 to 5 days

**Questions**

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.