

Lymph3Cx Assay, Primary Mediastinal Large B-cell Lymphoma and Diffuse Large B-cell Lymphoma, mRNA Gene Expression, NanoString, Tissue

Test ID: PM3CX

Explanation: The specimen requirements have been updated as shown below.

Current Specimen Required
<p>Submit only 1 of the following specimens:</p> <p>Preferred: Specimen Type: Tissue slides Slides: 1 stained and 7 unstained</p> <p>Collection Instructions: Submit 1 slide stained with hematoxylin and eosin and 7 consecutive, unstained, 5-micron thick sections placed on positively charged slides.</p> <p>Additional Information: Paraffin embedded specimens can be from any anatomic location (skin, soft tissue, lymph node, etc).</p> <p>Acceptable: Specimen Type: Tissue block Collection Instructions: Submit a formalin-fixed, paraffin-embedded tumor tissue block.</p>

New Specimen Required
<p>Submit only 1 of the following specimens:</p> <p>Preferred: Specimen Type: Tissue slides Slides: 1 Hematoxylin and eosin stained and 6 to 10 unstained</p> <p>Collection Instructions:</p> <ol style="list-style-type: none">1. Submit 1 slide stained with hematoxylin and eosin.2. If the formalin-fixed paraffin-embedded (FFPE) biopsy specimen is a core needle biopsy or fine needle aspirate, submit 10 unstained slides.3. If the FFPE biopsy specimen is a resection, submit 6 unstained slides.4. Submitted slides should be 4- to 5-micron thick sections placed on positively charged slides. <p>Additional Information: Paraffin embedded specimens can be from any anatomic location (skin, soft tissue, lymph node, etc).</p> <p>Acceptable: Specimen Type: Tissue block Collection Instructions: Submit an FFPE tumor tissue block.</p>

Questions

Contact Melissa Lonzo, Laboratory Resource Coordinator at 800-533-1710.