



Collagen Type II antibody

Test ID: FFTYC

Explanation: Per notification received from IMMCO Diagnostics, effective immediately, specimen stability for Test ID FFTYC will experience the following change.

Current Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	365 days
	Ambient	7 days

New Specimen Stability		
Specimen	Temperature	Time
Serum	Refrigerated (preferred)	7 days
	Frozen	365 days

Questions:

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.