



Amantadine, Serum or Plasma

Test ID: FAMAN

Explanation: Effective immediately, test ID FAMAN, referred to Medtox, will reflect the following changes:

Current Specimen Stability		
Specimen	Temperature	Time
Varies	Refrigerated (preferred)	14 days
	Ambient	72 hours
	Frozen	180 days

New Specimen Stability		
Specimen	Temperature	Time
Varies	Refrigerated (preferred)	21 days
	Ambient	21 days
	Frozen	28 days

Questions:

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.