

Fungitell, BAL

Test ID: FUNBL

Explanation: Effective immediately, test ID FUNBL, referred to Eurofins Viracor, will reflect the following change.

Current Specimen Stability		
Specimen	Temperature	Time
Lavage	Frozen (preferred)	30 days
	Ambient	4 days
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
Lavage	Frozen (preferred)	14 days
	Ambient	14 days
	Refrigerated	14 days

Questions:

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.