

Leukotriene E4, Random, Urine

Test ID: RLTE4

Explanation: The following stability changes will occur on the effective date.

Current Specimen Stability		
Specimen	Temperature	Time
Urine	Frozen (preferred)	28 days
	Ambient	24 hours
	Refrigerated	7 days

New Specimen Stability		
Specimen	Temperature	Time
Urine	Frozen (preferred)	28 days
	Ambient	8 hours
	Refrigerated	72 hours

Questions

Contact Joshua Yang, Laboratory Resource Coordinator at 800-533-1710.