

***Yersinia* Culture, Feces**

Test ID: YERSC; performed at Mayo Clinic Laboratories Florida.

Explanation:

Due to manufacturer supply shortage of culture media, YERSC will be placed test down effective immediately. Specimens received in the lab will be cancelled and a remedy will be submitted for the recommended alternative test below. The duration of the test down is expected to last 4 weeks or more. Another update will be provided once standard testing resumes.

Recommended Alternative Test:

Gastrointestinal Pathogen Panel, PCR, Feces

Test ID: GIP; performed at Mayo Clinic Laboratories Florida.

Methodology:

Multiplex Polymerase Chain Reaction (PCR)

Reference Values:

Negative (for all targets)

Specimen Requirements:

Container/Tube:

Preferred: Specific modified Cary-Blair transport system; see Additional Information for acceptable collection media

Acceptable: Approved Cary-Blair transport system (15 mL of non-nutritive transport medium containing phenol red as a pH indicator)

Specimen Volume: Representative portion of feces

Collection Instructions:

1. Collect fresh fecal specimen and submit 1 gram or 5 mL in container with transport medium.
2. Within 2 hours of collection, place feces in preservative.
3. Submit preserved feces in original container. **Do not aliquot.**

4. If unpreserved specimens are received, testing will be canceled.

Minimum Volume: 1 mL

Specimen Stability Information:

Specimen Type	Temperature	Time
Fecal	Ambient (preferred)	4 days
	Refrigerated	4 days

CPT Code:

87507

Day(s) Performed: Monday through Sunday **Report Available:** 1 to 2 days

Questions

Contact Betheney Wiles, Laboratory Resource Coordinator at 800-533-1710.