

Cytokine Panel 13, CSF

Test ID: FCP13

Explanation: Effective October 19, 2026, test ID FCP13, referred to ARUP Laboratories, will reflect the following changes:

Current Specimen Required	New Specimen Required
Container/Tube: Sterile vial Specimen Volume: 1 mL Collection Instructions: 1. Transfer 1 mL CSF to a sterile plastic container. 2. Send frozen. Note: Critical frozen. Additional specimens must be submitted when multiple tests are ordered.	Container/Tube: Sterile vial Specimen Volume: 1 mL Collection Instructions: 1. Transfer 1 mL CSF to a sterile plastic container. 2. Freeze within 30 minutes and send frozen. Note: Critical frozen. Additional specimens must be submitted when multiple tests are ordered.

Questions:

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.