

Reporting Title: HIV-1/-2 Ag and Ab Screen, P**Performing Location** Rochester**Ordering Guidance:**

If specimen is from either autopsy or cadaver blood sources, the proper FDA-licensed assay is HV1CD / HIV-1 and HIV-2 Antibodies for Cadaveric or Hemolyzed Specimens, Serum.

Screening, supplemental or confirmatory serologic tests for HIV-1 or HIV-2 antibodies cannot distinguish between active neonatal HIV infection and passive transfer of maternal HIV antibodies in infants during the postnatal period (up to 2 years old). Diagnosis of HIV infection in newborns and infants up to 2 years old should be made by virologic tests, such as detection of HIV-1 DNA and RNA (HIVP / HIV-1 DNA and RNA Qualitative Detection by PCR, Plasma) or HIV-1 RNA (HIVQN / HIV-1 RNA Detection and Quantification, Plasma).

New York State clients: This test should not be requested for maternal/newborn HIV screening on specimens originating in New York State, due to state regulatory requirements for expedited result reporting.

Specimen Requirements:**Collection Container/Tube:** Lavender top (EDTA)**Submission Container/Tube:** Plastic vial**Specimen Volume:** 4 mL**Collection Instructions:**

1. Centrifuge blood collection tube per collection tube manufacturer's instructions (eg, centrifuge and aliquot within 2 hours of collection for BD Vacutainer tubes).
2. Aliquot plasma into a plastic vial.

Forms:

If not ordering electronically, complete, print, and send a [General Request](#) (T239) with the specimen.

Specimen Type	Temperature	Time	Special Container
Plasma	Refrigerated (preferred)	6 days	
	Frozen	30 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC
HIVC3	HIV-1/-2 Ag and Ab Screen, P	Alphanumeric		56888-1

LOINC and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

87389

G0475

Reflex Tests:

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
HVDIP	HIV Ab Confirm / Differentiation, P	1	86701	No	Yes
HIVQN	HIV-1 RNA Detect/Quant, P	1	87536	No	Yes

Result Codes for Reflex Tests:

Test ID	Result ID	Reporting Name	Type	Unit	LOINC
HIVQN	113581	HIV-1 RNA Detect/Quant, P	Alphanumeric	copies/mL	70241-5
HVDIP	91947	HIV-1 Ab Differentiation, P	Alphanumeric		68961-2
HVDIP	91951	HIV-2 Ab Differentiation, P	Alphanumeric		81641-3

Reference Values:

Negative