

**Reporting Title:** Hemoquant, F**Performing Location** Rochester**Specimen Requirements:****Patient Preparation:** Patient should refrain from ingesting red meat and aspirin-containing products (eg, Excedrin, Aspirin) for 3 days prior to specimen collection.**Collection Container/Tube:** Spoon-like sampler from kit (T134)**Submission Container/Tube:** Screw-capped tube**Specimen Volume:** 1 g**Collection Instructions:** Collect random specimen from a single defecation.

Specimen Type	Temperature	Time	Special Container
Fecal	Refrigerated (preferred)	7 days	
	Ambient	7 days	
	Frozen	0 hours	

**Result Codes:**

Result ID	Reporting Name	Type	Unit	LOINC
2410	Fecal Hemoglobin	Numeric	mg Hb/g	27396-1

LOINC and CPT codes are provided by the performing laboratory.

**Supplemental Report:**

No

**CPT Code Information:**

84126

**Reference Values:**

Normal:

&lt; or =2.0 mg total hemoglobin/g feces

Marginal:

2.1-4.0 mg total hemoglobin/g feces\*

\*2.1-4.0 mg Hb/g is considered marginally elevated, but not clinically significant, if red meat, warfarin, or aspirin was ingested 72 hrs prior to collection.

Elevated:

&gt;4.0 mg total hemoglobin/g feces