

Reporting Title: Targeted Stimulant Screen, U**Performing Location** Rochester**Specimen Requirements:**

Only orderable as part of profile. For more information see:

-CSMPU / Controlled Substance Monitoring Panel, Random, Urine

-TSPU / Targeted Stimulant Screen, Random, Urine

Supplies: Aliquot Tube, 5 mL (T465)**Collection Container/Tube:** Plastic urine container**Submission Container/Tube:** Plastic, 5-mL vial**Specimen Volume:** 1 mL**Collection Instructions:**

1. Collect a random urine specimen.
2. No preservative

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	14 days	
	Frozen	28 days	

Supplemental Report:

No

CPT Code Information:

80326

G0480 (if appropriate)

Reference Values:

Only orderable as part of profile. For more information see:

-CSMPU / Controlled Substance Monitoring Panel, Random, Urine

-TSPU / Targeted Stimulant Screen, Random, Urine

Not detected

Cutoff concentrations:

Methamphetamine: 100 ng/mL

Amphetamine: 100 ng/mL

3,4-Methylenedioxymethamphetamine (MDMA): 100 ng/mL

3,4-Methylenedioxy-N-ethylamphetamine (MDEA): 100 ng/mL

3,4-Methylenedioxyamphetamine (MDA): 100 ng/mL

Ephedrine: 100 ng/mL

Pseudoephedrine: 100 ng/mL

Phentermine: 100 ng/mL

Phencyclidine (PCP): 20 ng/mL

Methylphenidate: 20 ng/mL

Ritalinic acid: 100 ng/mL