

**Reporting Title:** 11-Desoxycortisol**Performing Location:** Esoterix Endocrinology**Specimen Requirements:****Collection Container/Tube:** Red-Top**Submission Container/Tube:** Plastic vial**Specimen Volume:** 1 mL**Acceptable:** SST**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Serum must be separated from cells within 45 minutes of venipuncture. Spin down and send 1 mL of serum frozen in a plastic vial.

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	0 hours	
	Ambient	72 hours	
	Refrigerated	72 hours	

**Supplemental Report:**

No

**CPT Code Information:**

82634

**Reference Values:**

Age	Range (ng/dL)
Premature (26 - 28 weeks) Day 4	110-1376
Premature (31 - 35 weeks) Day 4	48-579
Newborn Day 3	13-147
Newborn 1 to 11 months	<10-156
Prepubertal 8 AM	20-155
Prepubertal Children and Adults 8 AM	12-158