

# **Test Definition: TTF40**

TTF41 (SPT24) + p40 Immunostain, Technical Component Only

Reporting Title: TTF1 (SPT24) + p40 IHC, Tech Only

Performing Location: Rochester

### **Ordering Guidance:**

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

## **Shipping Instructions:**

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

### **Specimen Requirements:**

Specimen Type: Tissue

**Supplies:** Immunostain Technical Only Envelope (T693) **Container/Tube:** Immunostain Technical Only Envelope

**Preferred:** 

-Formalin-fixed, paraffin-embedded tissue block

OR

-2 Unstained, positively charged glass slides (25- x 75- x 1-mm) per test ordered; sections 4-microns thick

Acceptable: None

#### Forms:

If not ordering electronically, complete, print, and send an <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains</u>
<u>Request</u> (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
602648	TTF1 (SPT24) + p40 IHC, Tech Only	Alphanumeric		Bill only; no result

LOINC® and CPT codes are provided by the performing laboratory.

### **Supplemental Report:**

No

#### **CPT Code Information:**

88344-TC

#### **Reflex Tests:**



**Test Definition: TTF40** 

TTF41 (SPT24) + p40 Immunostain, Technical Component Only

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
OTMHI	IHC Multiplex, Tech Only	1	88344	No	No

### **Reference Values:**

This is not an orderable test. Order PATHC / Pathology Consultation. The consultant will determine the need for special stains.