

Test Definition: SGUR

Specific Gravity, Random, Urine

Reporting Title: Specific Gravity, U

Performing Location: Mayo Clinic Laboratories - Rochester Main Campus

Ordering Guidance:

Urine with contrast dye, glucose, or excessive protein should not be evaluated with this test. In these cases, urine osmolality is a better measure of urine concentration. Order UOSMU / Osmolality, Random, Urine.

Specimen Requirements:

Container/Tube: Plastic urine container

Specimen Volume: 20 mL **Collection Instructions:**

1. Collect a random urine specimen.

2. No preservative.

Forms:

If not ordering electronically, complete, print, and send a Renal Diagnostics Test Request (T830) with the specimen.

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	7 days	
	Frozen	7 days	

Result Codes:

Result ID	Reporting Name	Туре	Unit	LOINC®
SGUR	Specific Gravity, U	Numeric		In Process

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

Reference Values:

1.002-1.030