

## **Test Definition: PCEAL**

Carcinoembryonic Antigen, Polyclonal Immunostain, Technical Component Only

Reporting Title: P-CEA IHC, Tech Only

Performing Location: Rochester

### Ordering Guidance:

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

## **Shipping Instructions:**

Attach the green "Attention Pathology" address label (T498and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

## **Specimen Requirements:**

Specimen Type: Tissue

**Supplies:** Immunostain Technical Only Envelope (T693) **Container/Tube:** Immunostain Technical Only Envelope

**Preferred:** 

-Formalin-fixed, paraffin-embedded tissue block

OR

-2 Unstained, positively charged glass slides (25- x 75- x 1-mm) per test ordered; sections 4-microns thick

Acceptable: None

#### Forms:

If not ordering electronically, complete, print, and send a <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains</u>
<u>Request</u> (T763) with the specimen.

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

## **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
70847	P-CEA IHC, Tech Only	Alphanumeric		Bill only; no result

LOINC® and CPT codes are provided by the performing laboratory.

## **Supplemental Report:**

No

### **CPT Code Information:**

88342-TC, Primary 88341-TC, If additional IHC



# **Test Definition: PCEAI**

Carcinoembryonic Antigen, Polyclonal Immunostain, Technical Component Only

### **Reflex Tests:**

Test Id	Reporting Name	CPT Units	CPT Code	Always Performed	Available Separately
IHTOI	IHC Initial, Tech Only	1	88342	No	No
IHTOA	IHC Additional, Tech Only	1	88341	No	No