

## **Test Definition: PBNA**

Lead, Nails

**Reporting Title:** Lead, Nails **Performing Location:** Rochester

#### **Necessary Information:**

Indicate source of nails (fingernails or toenails), if known

#### **Specimen Requirements:**

Supplies: Hair and Nails Collection Kit (T565)

**Source**: Fingernails or toenails **Specimen Volume**: 0.2 g **Collection Instructions**:

1. Prepare and transport specimen per the instructions in the kit or see Collecting Hair and Nails for Metals Testing.

2. Clippings should be taken from all 10 fingernails or toenails.

Specimen Type	Temperature	Time	Special Container
Nail	Ambient (preferred)		
	Refrigerated		
	Frozen		

### Ask at Order Entry (AOE) Questions:

Test ID	Question ID	Description	Туре	Reportable
PBNA	PBNSC	Specimen Source:	Answer List	Yes
		• Finger		
		• Toe		
		Finger and Toe		

### **Result Codes:**

Result ID	Reporting Name	Туре	Unit	LOINC®
2506	Lead, Nails	Numeric	mcg/g	8202-4
PBNSC	Specimen Source	Alphanumeric		31208-2

LOINC® and CPT codes are provided by the performing laboratory.

#### **Supplemental Report:**

No

#### **CPT Code Information:**

83655

### **Reference Values:**

<4.0 mcg/g of nails



# **Test Definition: PBNA**

Lead, Nails

Reference values apply to all ages.