


HIV Testing Algorithm (Fourth-Generation Screening Assay), Including Follow-up of Reactive Rapid Serologic Test Results



* Plasma or serum specimens from individuals with reactive results from rapid HIV antigen/antibody, 3rd generation HIV-1/-2 antibody, or 5th generation HIV antigen/antibody screening tests can enter this testing algorithm at [HVDIP](#) or [HIVDI](#) with applicable downstream testing.

¹ Submit new plasma specimen for [HIVQN](#).

² Contact state department of health HIV surveillance unit for arrangement to send additional specimens to CDC for HIV-2 testing.