

## **Test Definition: FLNDR**

Flounder (Bothidae/Pleuronectidae Fam) IgE

## **Overview**

#### **Method Name**

ImmunoCAP FEIA

#### **NY State Available**

Yes

## **Specimen**

## **Specimen Type**

Serum

## **Specimen Required**

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

## **Specimen Minimum Volume**

0.5 mL

## **Reject Due To**

Hemolysis	Mild OK; Gross reject
Lipemia	Mild OK; Gross reject
Icterus	Mild OK; Gross reject
Other	NA

## **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

## **Clinical & Interpretive**

## **Reference Values**

<0.35 kU/L



## **Test Definition: FLNDR**

Flounder (Bothidae/Pleuronectidae Fam) IgE

#### Interpretation

<u>Class</u>	<u>lgE (kU/L)</u>	<u>Comment</u>
0	<0.10	Negative
0/1	0.10 – 0.34	Equivocal/Borderline
1	0.35 – 0.69	Low Positive
2	0.70 – 3.49	<b>Moderate Positive</b>
3	3.50 – 17.49	High Positive
4	17.50 – 49.99	Very High Positive
5	50.00 – 99.99	Very High Positive
6	>99.99	Very High Positive

#### **Performance**

## **PDF Report**

No

## Day(s) Performed

Monday through Friday

#### **Report Available**

1 to 6 days

## **Performing Laboratory Location**

**Eurofins Viracor** 

## **Fees & Codes**

#### **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

#### **Test Classification**

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

#### **CPT Code Information**

86003

#### **LOINC®** Information

Test ID	Test Order Name	Order LOINC® Value



# **Test Definition: FLNDR**

Flounder (Bothidae/Pleuronectidae Fam) IgE

FLNDR	Flounder IgE	7316-3
Result ID	Test Result Name	Result LOINC® Value
Z3815	Flounder IgE	7316-3
Z3816	CLASS	102311-8