

Test Definition: ZW191

Miscellaneous University of Texas Health
Center at Tyler Microbiology

Overview

Special Instructions

• Infectious Specimen Shipping Guidelines

Method Name

Varies

NY State Available

No

Specimen

Specimen Type

Varies

Specimen Required

Varies

This is a Miscellaneous Referral Test. For specific requirements, please refer to the referred tests list under the catalog tab on the MCL Website (external) or the Referral Catalog (CRM) (internal). If unable to find the correct test, request forms, or specimen requirements, contact Customer Service at 1-800-533-1710 or 6-5700.

NOTE: Provide when ordering

- 1. Test name
- 2. Performing lab code
- 3. Specimen Type

Reject Due To

Specimens	Varies
other than	
Anticoagulants	NA NA
other than	
Hemolysis	NA
Lipemia	NA
Icteric	NA

Specimen Stability Information



Test Definition: ZW191

Miscellaneous University of Texas Health
Center at Tyler Microbiology

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

Clinical & Interpretive

Reference Values

Test Performed By: University of Texas Health Center

Department of Microbiology

11937 US Hwy 271 Tyler, TX 75708

Performance

PDF Report

Referral

Day(s) Performed

Varies

Report Available

Varies

Performing Laboratory Location

University of Texas Health Center

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

CPT Code Information

Varies

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
ZW191	Misc Univ TX-Tyler Micro	51991-8



Test Definition: ZW191

Miscellaneous University of Texas Health Center at Tyler Microbiology

Result ID	Test Result Name	Result LOINC® Value
ZT191	Test Name	19145-2
ZR191	Result	19146-0
ZF191	Flag	No LOINC Needed
ZV191	Reference Value	19147-8
ZU191	Unit of Measure	No LOINC Needed