

Protein, Total, 24 Hour, Urine

Test ID: PTU

Explanation: On the effective date, reporting name, specimen stability and result codes will be updated as indicated below.

Current Reporting Name
Protein, Total, U

New Reporting Name
Protein, Total, 24 HR, U

Current Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerate (preferred)	14 days
	Ambient	24 hours
	Frozen	14 days

New Specimen Stability		
Specimen	Temperature	Time
Urine	Refrigerate (preferred)	14 days
	Ambient	24 hours
	Frozen	30 days

Current Result Codes	
Result ID	Reporting Name
TP2	Total Protein
TM23	Collection Duration
VL21	Urine Volume
CONC1	Concentration

New Result Codes	
Result ID	Reporting Name
TP2	Total Protein, 24 HR, U
TM23	Collection Duration
VL21	Urine Volume

Questions

Contact Nancy Benson, Laboratory Technologist Resource Coordinator at 800-533-1710.