
Diuretic Profile, Urine

Explanation: The following new orderable, referred to Medtox Laboratory, will be available on August 25, 2025.

Test ID: FDISU

Methods:

High Performance Liquid Chromatography with Ultraviolet Detection; (HPLC-UV)

Specimen Requirements:

Specimen Type: Urine

Container/Tube: Plastic preservative-free urine container

Specimen Volume: 10 mL

Collection Instructions:

1. Collect 10 mL random urine without preservative.
2. Send specimen refrigerated

Specimen Stability Information:

Specimen Type	Temperature	Time
Urine	Refrigerated (preferred)	14 days
	Ambient	72 hours
	Frozen	14 days

Fee:
\$130.83

CPT Code:
80377

Day(s) Setup: Thursday

Report Available: 7 - 9 day

Questions

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.