

Glutamic Acid Decarboxylase (GAD) Autoantibody

Explanation: The following new orderable, referred to LabCorp Burlington, will be available on October 13, 2026.

Test ID: FGADA

Method:

Enzyme-Linked Immunosorbent Assay (ELISA)

Reference Values:

0.0-5.0 units/mL

Specimen Requirements:

Collection Container/Tube:

Preferred: Red top

Acceptable: Serum gel

Submission Container/Tube: Plastic vial

Specimen Volume: 0.5 mL Serum

Collection Instructions:

1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Send refrigerated.

Specimen Minimum Volume

Serum: 0.4 mL

Note: This volume does not allow for repeat testing.

Specimen Stability Information:

Specimen Type	Temperature	Time
Serum	Refrigerated (Preferred)	14 days
	Ambient	14 days
	Frozen	14 days

CPT Code:

86341

Day(s) Setup: Monday through Saturday

Report Available: 5 to 8 days

Questions:

Contact MCS Referrals Supervisors Amy Bluhm or Heather Hilgart at 800-533-1710.