Test Definition: POSA

Posaconazole, Serum

Reporting Title: Posaconazole, S
Performing Location: Mayo Clinic Laboratories - Rochester Superior Drive

Specimen Requirements:

Supplies: Sarstedt Aliquot Tube, 5 mL (T914)
Collection Container/Tube: Red top (serum gel/SST are **not acceptable**)
Submission Container/Tube: Plastic vial
Specimen Volume: 2 mL serum
Collection Instructions: Within 2 hours of collection, centrifuge, and aliquot serum into plastic vial.

Forms:

If not ordering electronically, complete, print, and send a Therapeutics Test Request (T831) with the specimen.

Specimen Type	Temperature	Time	Special Container
Serum Red	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Result Codes:

Result ID	Reporting Name	Type	Unit	LOINC®
89591	Posaconazole, S	Numeric	ng/mL	53731-6

LOINC® and CPT codes are provided by the performing laboratory.

Supplemental Report:

No

CPT Code Information:

80187

Reference Values:

>700 ng/mL (trough)